



10000 SERIES

POWDER FREE LATEX EXAMINATION GLOVES – TEXTURED

Our most popular powder free exam glove provides exceptional comfort and quality at a competitive price point without sacrificing consistent dependability. This glove has a micro-textured grip and a precise fit providing reduced hand fatigue.



GLOVE SPECIFICATIONS

| WEIGHT & DIMENSION | | | |
|--------------------|---------|-------------|------------|
| Reorder# | Size | Length | Palm Width |
| 10110 | X-Small | 238 - 242mm | 75 ± 5mm |
| 10111 | Small | 238 - 242mm | 85 ± 5mm |
| 10112 | Medium | 238 - 242mm | 95 ± 5mm |
| 10113 | Large | 238 - 242mm | 105 ± 5mm |
| 10114 | X-Large | 238 - 242mm | 115 ± 5mm |

| MIL THICKNESS (SINGLE WALL) | | |
|-----------------------------|---------|---------|
| | Palm | Finger |
| All Sizes | 4.0 mil | 4.5 mil |

| PHYSICAL PROPERTIES | |
|---------------------------------|--------------|
| | Before Aging |
| Tensile Strength Mpa (minimum) | 18 |
| Ultimate Elongation % (minimum) | 650 |

| PERFORMANCE REQUIREMENTS | | | |
|--------------------------|--|----------------------------------|-----|
| Characteristic | Related Defect | Inspection Level | AQL |
| Freedom from holes | Holes, lumps | G - 1 / S - 4 Multiple Normal | 1.5 |
| Physical dimension | Length, width and thickness | S - 2 | 4.0 |
| Physical properties | Before aging and after accelerated aging | S - 2 | 4.0 |
| Visual defect | Stains etc. | S - 4 | 2.5 |

| PACKAGING | |
|---------------|----------------------------|
| 10 boxes/case | 100 gloves/box (by weight) |