

KIMTECH™ Purple Nitrile™ Ambidextrous Exam Glove, 9.5”

Product Information

Code	Size	Palm Width*(mm)	Length* (mm)	Quantity
55080	XS	70	242	100 gloves/box; 10 boxes/case
55081	S	80	242	100 gloves/box; 10 boxes/case
55082	M	95	242	100 gloves/box; 10 boxes/case
55083	L	110	242	100 gloves/box; 10 boxes/case
55084	XL	120	242	90 gloves/box; 10 boxes/case

*As measured by ASTM D 6319; 2.5 AQL S-2¹



Property	Description
General Product Information	Cleared for use in chemotherapy; tested for use with 29 chemotherapy drugs per ASTM D6978-05
General Description	Purple nitrile, powder-free, exam, non-sterile, textured fingertips, ambidextrous gloves, beaded cuff, and single chlorination, 9.5". Not made with natural rubber latex.

Physical Property	Unit	Upper Limit	Test Method
Watertightness	NA	Pass	ASTM D 5151 (1.0 AQL)
Residual Powder	mg/glove	2.0	ASTM D 6124 (2.5 AQL)

Thickness	Target	Target	Test Method (2.5 AQL S-2)
Middle Finger	5.91 mil	0.15 mm	ASTM D 6319
Palm	4.72 mil	0.12 mm	ASTM D 6319
Cuff	3.54 mil	0.09 mm	ASTM D 6319

Property	Tensile Strength	Ultimate Elongation	Median Force at Break
Before Aging	21.0 MPa	550%	11 Newtons
After Aging	21.0 MPa	500%	9 Newtons
Testing Method (2.5 AQL S-2)	ASTM D 412 ASTM D 6319 ASTM D 573	ASTM D 412 ASTM D 573 ASTM D 6319	EN 455-2

¹ AQL levels defined according to ISO 2859-1 (current edition) Sampling Procedures for Inspection by Attributes.

Workmanship Requirement:

Gloves shall be made in accordance with reasonable industry practice with respect to defects, dirt and contamination. Gloves should be packaged as specified having described size and color in the appropriate quantity as listed. Shipping Case and inner packaging shall be labeled with product name, size, lot/batch and catalog number.

For Additional information: Contact your Kimberly-Clark representative, our Customer Service Team (888)-346-GOKC (4652) or KCPIInfo@kcc.com

Rev. February 24, 2022

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products. **THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

© Registered Trademark or * Trademark of Kimberly-Clark Worldwide, Inc. Marques déposées de Kimberly-Clark Worldwide, Inc. ou de ses filiales. © KCWW.

